

EXTENDS *Naturals*

VARIABLES *count, pc, target, photon*

CONSTANT *Intensity, Time, Mode*

X ray intensity is set to a constant 25 *MeV*

$XRayIntensity \triangleq 25$

In “Electron” mode the electron beam intensity is lesser than the current in “XRay” mode.

$TypeInvariant \triangleq \wedge Time \in Nat \wedge Time \neq 0$
 $\wedge Mode \in \{“Electron”, “XRay”, “Done”\}$
 $\wedge Intensity \leq XRayIntensity$

The “count” variable is decreased at each step and it stops when it is zeroed. In the “XRay” mode, the target must be placed in the path of the beam to act as a filter and collimator. In the “Electron” mode, the target is not needed due to the low ‘lightning’ intensity.

$Init \triangleq \wedge count = Time$
 $\wedge pc = Mode$
 $\wedge IF pc = “XRay” THEN \wedge photon = XRayIntensity$
 $\wedge target = “On”$
 $ELSE \wedge photon = Intensity$
 $\wedge target = “Off”$

Electron mode

$Treat1 \triangleq \wedge pc = “Electron”$
 $\wedge target = “Off”$
 $\wedge IF count = 0 THEN \wedge pc' = “Done”$
 $\wedge photon' = 0$
 $\wedge UNCHANGED \langle target \rangle$
 $ELSE \wedge count' = count - 1$
 $\wedge UNCHANGED \langle pc, photon, target \rangle$

XRay mode

$Treat2 \triangleq \wedge pc = “XRay”$
 $\wedge target = “On”$ Confirm that the target is activated.
 $\wedge IF count = 0 THEN \wedge pc' = “Done”$
 $\wedge photon' = 0$
 $\wedge target' = “Off”$
 $ELSE \wedge count' = count - 1$
 $\wedge UNCHANGED \langle pc, photon, target \rangle$

$Next \triangleq Treat1 \vee Treat2$

$vars \triangleq \langle count, pc, target, photon \rangle$

$Spec \triangleq Init \wedge \Box [Next]_{\langle vars \rangle}$